



## Co-Dx Logix Smart™ SARS-CoV-2 (genes RdRp/E)

Innovating Molecular Diagnostic Solutions

The Logix Smart™ SARS-CoV-2 (genes RdRp/E) test is a real-time RT-PCR multiplex test intended for the in vitro qualitative detection of nucleic acid from severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), targeting the genes RdRp in the polygene Orf1ab region and gene E of the virus genome, in lower respiratory tract samples (e.g., bronchoalveolar lavage, sputum, tracheal aspirate), upper respiratory tract samples (e.g., nasopharyngeal and oropharyngeal swabs), and saliva from individuals suspected of COVID-19, as determined by their healthcare provider.

Results are used for the identification of SARS-CoV-2 RNA during the acute phase of infection. The Logix Smart™ SARS-CoV-2 (genes RdRp/E) is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and in vitro diagnostic procedures. The test kit has been tested with the QIAamp Viral RNA Mini Kit (Qiagen), Sbeadex Viral RNA Purification kit (Biosearch Technologies), Viral DNA/RNA kit (CW Bio), and HighPrep Viral DNA/RNA kit (MagBio) extraction systems on the Co-Dx Box™ (Co-Diagnostics, Inc.) Mic qPCR Cyclers (BMS, Biomolecular Systems), QuantStudio 5 (Thermo Fisher Scientific), CFX96 (Bio-Rad).

*For in vitro diagnostic use. For professional use only.*

### Co-Dx Logix Smart™ SARS-CoV-2

- Regulatory Status: European Union (IVD)
- Includes internal control to verify sample quality
- Includes a positive control to verify master mix quality
- Produces results that are easy to interpret
- For use with lower and upper respiratory tract specimen as well as saliva samples



<b>Intended Use</b>	Qualitative real-time RT-PCR multiplex test for in vitro detection of nucleic acid from severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2).
<b>Sample Type</b>	Lower respiratory samples (e.g., bronchoalveolar lavage, sputum, tracheal aspirate), upper respiratory samples (e.g., nasopharyngeal and oropharyngeal swabs), and saliva.
<b>User</b>	Qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and in vitro diagnostic procedures.
<b>Limit of Detection</b>	The Limit of Detection (LoD) for Logix Smart™ SARS-CoV-2 utilizing QIAamp RNA Viral Mini Kit (cat. no. 52904/ 52906, Qiagen) was confirmed to be 0.8 copies/μL. For the HighPrep Viral DNA/RNA Kit (MagBio, CAT#HPV-DR96) the LoD was confirmed to be 1.0 copies/μL. For the Viral DNA/RNA Kit (CW Bio, CAT#CW3126M) the LoD was confirmed to be 2.0 copies/μL. For the sbeadex Viral RNA Purification Kit was run on the oKtopure (Biosearch Technologies, CAT#NAP-40-026-04) the LoD was confirmed to be 6.0 copies/μL.
<b>Sensitivity*</b>	99.448%
<b>Specificity*</b>	98.908%
<b>Clinical Matrix used for analytical verification</b>	Lower respiratory samples (e.g., bronchoalveolar lavage, sputum, tracheal aspirate), upper respiratory samples (e.g., nasopharyngeal and oropharyngeal swabs), and saliva.
<b>Analytical Specificity (in silico analysis)</b>	No microorganism in the in silico analysis has revealed significant homology between the cross-reactivity microorganisms, including the ones of relevance listed below.  IT DOES NOT cross-react with the following microorganisms: Human coronavirus 229E, Adenovirus, Influenza C, Human coronavirus OC43, Human Metapneumovirus (hMPV), Parechovirus, Human coronavirus HKU1, Parainfluenza virus 1-4, Corynebacterium diphtheriae, Human coronavirus NL63, Influenza A & B, Legionella non-pneumophila, SARS-coronavirus, Enterovirus, Bacillus anthracis (Anthrax), MERS- coronavirus, Respiratory syncytial virus, Moraxella catarrhalis, Rhinovirus, Neisseria elongata, Chlamydia pneumoniae, Neisseria meningitides, Haemophilus Influenza, Leptospiriosis, Legionella pneumophila, Chlamydia psittaci, Mycobacterium tuberculosis, Coxiella burnetii (Q-Fever), Streptococcus pneumoniae, Staphylococcus aureus, Streptococcus pyogenes, Bordetella pertussis, Mycoplasma pneumoniae, Pneumocystis jirovecii (PJP), Pooled human nasal wash – to represent diverse microbial flora in the human respiratory tract, Candida albicans, Pseudomonas aeruginosa, Staphylococcus epidermidis, and Staphylococcus salivarius.
<b>Time to detection</b>	Approximately 90 minutes, depending on the instrument used.
<b>Extraction System</b>	<ul style="list-style-type: none"> <li>• QIAamp Viral RNA Mini Kit (Qiagen)</li> <li>• Sbeadex Viral RNA Purification kit (Biosearch Technologies)</li> <li>• Viral DNA/RNA kit (CW Bio)</li> <li>• HighPrep Viral DNA/RNA kit (MagBio)</li> </ul>
<b>Thermal cycler compatibility</b>	<ul style="list-style-type: none"> <li>• Co-Dx Box™ Cycle (Co-Diagnostics, Inc.)</li> <li>• Mic qPCR Cycler (BMS, Biomolecular Systems)</li> <li>• QuantStudio 5 (Thermo Fisher Scientific)</li> <li>• CFX96 (Bio-Rad)</li> </ul> Each thermal cycler validated with the test and listed above is compatible with the following detection channels: <ul style="list-style-type: none"> <li>• FAM</li> <li>• CF560 (VIC)</li> <li>• CF610 (ROX)</li> </ul>

\* Sensitivity based on clinical study of 30 clinical remnant positive samples, 15 contrived Influenza A, and 15 contrived Influenza B, and 30 negative clinical samples.

## Each Co-Dx Logix Smart™ SARS-CoV-2 kit includes:

Cap Color	Component	Description	Amount
<b>Brown</b>	Co-Dx Logix Smart™ SARS-CoV-2 Master Mix	Proprietary blend of SARS-CoV-2 Co-Primers™ and PCR reagents	1×500 μL (100 reactions) or 1×1250 μL (250 reactions) or 1×25000 μL (5,000 reactions)
<b>Red</b>	Co-Dx Logix Smart™ SARS-CoV-2 Positive Control	Proprietary blend of SARS-CoV-2 synthetic templates	1×500 μL (100 reactions) or 1×1250 μL (250 reactions) or 1×25000 μL (5,000 reactions)
<b>Clear</b>	Nuclease Free Water	Water free of DNase/RNase activity	1×500 μL (100 reactions) or 1×1250 μL (250 reactions) or 1×25000 μL (5,000 reactions)

## Ordering Information:

Product Name	Product ID
Co-Dx Box™ Thermocycler	Request quote
Co-Dx Logix Smart™ Coronavirus Disease 2019 (COVID-19) CE-IVD Kit	COVID-K-002

